
Product Datasheet

COL17A1 Antibody (orb1262376)

Description

COL17A1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

This COL17A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 475-504 amino acids from the Central region of human COL17A1.

Target

COL17A1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25 For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


150 kDa

Uniprot ID
[Q9UMD9](#)
NCBI
[Q9UMD9](#)
Dilution Range


For FACS starting dilution is: 1:25 For WB starting dilution is: 1:2000

Expiration Date

12 months from date of receipt.



Overlay histogram showing A431 cells sta...



Western Blot at 1:2000 dilution Lane 1: ...