


Product Datasheet

TUBB1 Antibody (orb1262410)


Description	TUBB1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	This TBB1 antibody is generated from rabbits immunized with human TBB1 recombinant protein.
Target	TUBB1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q9H4B7
NCBI	Q9H4B7
Dilution Range	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot at 1:2000 dilution Lane 1: ...



Formalin-fixed and paraffin-embedded hum...



Flow cytometry analysis of MCF-7 cells (...)