
Product Datasheet

ADAMTS18 Antibody (orb1262459)

Description

ADAMTS18 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, IHC-P, WB

Immunogen

This ADAMTS18 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1066-1095 amino acids from the C-terminal region of human ADAMTS18.

Target

ADAMTS18

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



135 kDa

Uniprot ID
[Q8TE60](#)
NCBI
[Q8TE60](#)
Dilution Range


For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

Western blot analysis of ADAMTS18 Antibo...



Formalin-fixed and paraffin-embedded hum...



Western blot analysis in MDA-MB231 cell ...