


Product Datasheet

STK24 Antibody (orb1262858)

Description	STK24 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This MST3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 345-374 amino acids from the C-terminal region of human MST3.
Target	STK24
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q9Y6E0
NCBI	Q9Y6E0
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of MST3 Antibody i...



Formalin-fixed and paraffin-embedded hum...