


Product Datasheet

ALOX5AP Antibody (orb1262871)

Description	ALOX5AP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This ALOX5AP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 86-117 amino acids from the Central region of human ALOX5AP.
Target	ALOX5AP
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	P20292
NCBI	P20292
Dilution Range	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-ALOX5AP An...



Formalin-fixed and paraffin-embedded hum...