
Product Datasheet

MERTK Antibody (orb1262896)

Description

MERTK Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

This MERTK antibody is generated from rabbits immunized with a his tag recombinant protein of human MERTK.

Target

MERTK

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

 For WB starting dilution is: 1:1000
 For IF starting dilution is: 1:10~50
 For IHC-P starting dilution is: 1:50~10
 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


110 kDa

Uniprot ID
[Q12866](#)
NCBI
[Q12866](#)
Dilution Range


 For WB starting dilution is: 1:1000
 For IF starting dilution is: 1:10~50
 For IHC-P starting dilution is: 1:50~10
 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in Jurkat cell lin...



Fluorescent image of U251 cell stained w...



Formalin-fixed and paraffin-embedded hum...