
Product Datasheet

NTRK3 Antibody (orb1262998)

Description

NTRK3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This TrkC antibody is generated from rabbits immunized with a his tag recombinant protein of human TrkC.

Target

NTRK3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:100 For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


94 kDa

Uniprot ID
[Q16288](#)
NCBI
[Q16288](#)
Dilution Range


For IHC-P starting dilution is: 1:100 For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis in U87-MG cell lin...