




Product Datasheet

GPT Antibody (orb1263198)

Description	GPT Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This GPT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 86-112 amino acids from the N-terminal region of human GPT.	Western blot analysis in HL60 cell line ...
Target	GPT	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Formalin-fixed and paraffin-embedded hum...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50 For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	55 kDa	Flow cytometric analysis of 293 cells us...
Uniprot ID	P24298	
NCBI	P24298	
Dilution Range	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50 For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	