


Product Datasheet

VTN Antibody (orb1263201)


Description	VTN Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	This VTN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 352-379 amino acids from the C-terminal region of human VTN.
Target	VTN
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	P04004
NCBI	P04004
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of VTN Antibody in...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of NCI-H292 cel...