




Product Datasheet

GGT5 Antibody (orb1263256)

Description	GGT5 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This GGTLA1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 130-156 amino acids from the N-terminal region of human GGTLA1.	Western blot analysis of GGTLA1 antibody...
Target	GGT5	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	Western blot analysis of GGTLA1 antibody...
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	Formalin-fixed and paraffin-embedded hum...
Clonality	Polyclonal	
MW	62 kDa	
Uniprot ID	P36269	
NCBI	P36269	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	