
Product Datasheet

GALK2 Antibody (orb1263483)

Description

GALK2 Antibody


Species/Host

Rabbit

Reactivity


Human

Conjugation

Unconjugated

Tested

WB

Applications


Antibody is used in Western blot to dete...

Immunogen

This GALK2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 404-434 amino acids from the C-terminal region of human GALK2.

Target

GALK2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

50 kDa

Uniprot ID
[Q01415](#)
NCBI
[Q01415](#)
Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.