


Product Datasheet

GNAS Antibody (orb1263615)

Description	GNAS Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This GNAS antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 640-667 amino acids from the Central region of human GNAS.
Target	GNAS
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:4000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	111 kDa
Uniprot ID	Q5JWF2
NCBI	Q5JWF2
Dilution Range	For WB starting dilution is: 1:4000
Expiration Date	12 months from date of receipt.



Western Blot at 1:4000 dilution + Ramos ...