



Product Datasheet

DLGAP1 Antibody (orb1263626)

Description	DLGAP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This DLGAP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 162-190 amino acids from the N-terminal region of human DLGAP1.
Target	DLGAP1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	109 kDa
Uniprot ID	O14490
NCBI	O14490
Dilution Range	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.

Western blot analysis of DLGAP1 Antibody...



Flow cytometric analysis of MDA-231 cell...