


Product Datasheet

ARPC3 Antibody (orb1263766)

Description	ARPC3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This ARPC3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 131-158 amino acids from the C-terminal region of human ARPC3.
Target	ARPC3
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	O15145
NCBI	O15145
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of ARPC3 antibody ...



Formalin-fixed and paraffin-embedded hum...