
Product Datasheet

NPEPPS Antibody (orb1263795)

Description

NPEPPS Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This PSA antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 709-737 amino acids from the C-terminal region of human PSA.

Target

NPEPPS

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


103 kDa

Uniprot ID
[P55786](#)
NCBI
[P55786](#)
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Western blot analysis of PSA antibody in...



PSA Antibody IHC analysis in formalin fi...