


Product Datasheet

IL1RN Antibody (orb1263837)


Description	IL1RN Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	This IL1RN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 147-177 amino acids from the C-terminal region of human IL1RN.
Target	IL1RN
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P18510
NCBI	P18510
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot at 1:4000 dilution + 293 wh...



Western blot analysis in 293, HeLa cell ...



Flow cytometric analysis of HepG2 cells ...