



Product Datasheet

NRG1 Antibody (orb1263934)

Description	NRG1 Antibody	
Species/Host	Rabbit	Western blot analysis of anti-NRG1 Antib...
Reactivity	Human	
Conjugation	Unconjugated	 <p>Flow cytometric analysis of NCI-H460 cel...</p>
Tested Applications	FC, IHC-P, WB	
Immunogen	This NRG1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 198-229 amino acids from the Central region of human NRG1.	 <p>Formalin-fixed and paraffin-embedded hum...</p>
Target	NRG1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	70 kDa	
Uniprot ID	Q02297	
NCBI	Q02297	
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	