


Product Datasheet

RNF114 Antibody (orb1264104)

Description	RNF114 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western blot analysis in A2058 cell line...
Immunogen	This ZNF313 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 167-193 amino acids from the C-terminal region of human ZNF313.	
Target	RNF114	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	ZNF313 Antibody immunohistochemistry ana...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	26 kDa	
Uniprot ID	Q9Y508	
NCBI	Q9Y508	
Dilution Range	For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	