

Product Datasheet

TAF1D Antibody (orb1264129)

Description

TAF1D Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested
Applications

FC, IHC-P, WB

Immunogen

This TAF1D antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 62-90 amino acids from the N-terminal region of human TAF1D.

Target

TAF1D

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


32 kDa

Uniprot ID
Q9H5J8
NCBI
Q9H5J8
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in NCI-H292, Jurka...



TAF1D Antibody
immunohistochemistry
anal...



Flow cytometric
analysis of Jurkat
cells...