



Product Datasheet

PMCH Antibody (orb1264150)

Description	PMCH Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	Western blot analysis in K562 cell line ...
Immunogen	This PMCH antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 94-122 amino acids from the Central region of human PMCH.	
Target	PMCH	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	Western blot analysis in mouse brain tis...
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	PMCH Antibody immunohistochemistry analy...
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	19 kDa	
Uniprot ID	P20382	
NCBI	P20382	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	