
Product Datasheet

CTTNBP2 Antibody (orb1264163)

Description

CTTNBP2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This CTTB2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 72-100 amino acids from the N-terminal region of human CTTB2.

Target

CTTNBP2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


181 kDa


Uniprot ID
[Q8WZ74](#)
NCBI
[Q8WZ74](#)
Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



 Western blot analysis
in NCI-H460 cell I...

 CTTB2 Antibody
immunohistochemistry
anal...