


Product Datasheet

DHRS2 Antibody (orb1264223)

Description	DHRS2 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western Blot at 1:16000 dilution + HepG2...
Immunogen	This DHRS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 197-225 amino acids from the C-terminal region of human DHRS2.	
Target	DHRS2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	DHRS2 Antibody immunohistochemistry anal...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:16000 For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	30 kDa	
Uniprot ID	Q13268	
NCBI	Q13268	
Dilution Range	For WB starting dilution is: 1:16000 For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	