

Product Datasheet

RSAD1 Antibody (orb1264250)

Description

RSAD1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This RSAD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 299-328 amino acids from the C-terminal region of human RSAD1.

Target

RSAD1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


49 kDa

Uniprot ID
[Q9HA92](#)
NCBI
[Q9HA92](#)
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of RSAD1 Antibody ...



RSAD1 Antibody IHC analysis in formalin ...



Flow cytometric analysis of HepG2 cells ...