



Product Datasheet

ENPP7 Antibody (orb1264393)

Description	ENPP7 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This ENPP7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 286-313 amino acids from the Central region of human ENPP7.	Western blot analysis of ENPP7 Antibody ...
Target	ENPP7	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Flow cytometric analysis of HL-60 cells ...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	51 kDa	
Uniprot ID	Q6UWV6	
NCBI	Q6UWV6	
Dilution Range	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	