
Product Datasheet

FA2H Antibody (orb1264440)

Description	FA2H Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	Western blot analysis of FA2H Antibody i...
Tested Applications	IHC-P, WB	
Immunogen	This FA2H antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 121-148 amino acids from the Central region of human FA2H.	
Target	FA2H	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	Western blot analysis of FA2H Antibody i...
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	43 kDa	FA2H Antibody IHC analysis in formalin f...
Uniprot ID	Q7L5A8	
NCBI	Q7L5A8	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	