
Product Datasheet

APRT Antibody (orb1264524)

Description

APRT Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Immunogen

This APRT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 143-170 amino acids from the C-terminal region of human APRT.

Target

APRT

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


20 kDa

Uniprot ID
[P07741](#)
NCBI
[P07741](#)
Dilution Range


For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western Blot at 1:1000 dilution Lane 1: ...



Western blot analysis in MCF-7 cell line...



Formalin-fixed and paraffin-embedded hum...