

## Product Datasheet

### CAPZB Antibody (orb1264528)

## Description

CAPZB Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This CAPZB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 229-257 amino acids from the C-terminal region of human CAPZB.

### Target

CAPZB

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

31 kDa

### Uniprot ID

[P47756](#)

### NCBI

[P47756](#)


### Dilution Range

For WB starting dilution is: 1:1000


### Expiration Date

12 months from date of receipt.

m.brain  
72



Western blot analysis of CAPZB Antibody ...



Western blot analysis of CAPZB using rab...