


Product Datasheet

POLD1 Antibody (orb1264545)

Description	POLD1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This POLD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 710-739 amino acids from the Central region of human POLD1.
Target	POLD1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	P28340
NCBI	P28340
Dilution Range	For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of POLD1 Antibody ...



Flow cytometry analysis of Jurkat cells (...)