

Product Datasheet

U2AF1 Antibody (orb1264592)

Description

U2AF1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This U2AF1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 55-83 amino acids from the Central region of human U2AF1.

Target

U2AF1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


28 kDa

Uniprot ID
Q01081
NCBI
Q01081
Dilution Range


For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis of anti-U2AF1
Anti...



Western blot analysis of U2AF1 using rab...