
Product Datasheet

GADD45A Antibody (orb1264605)

Description

GADD45A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

This GADD45A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 133-165 amino acids from the C-terminal region of human GADD45A.

Target

GADD45A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



18 kDa

Uniprot ID
[P24522](#)
NCBI
[P24522](#)
Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of GADD45A Antibod...



Flow cytometric analysis of NCI-H460 cel...



Western blot analysis of GADD45A using r...