
Product Datasheet

AGL Antibody (orb1264727)

Description

AGL Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, WB

Applications
Immunogen

This AGL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1479-1510 amino acids from the C-terminal region of human AGL.

Target

AGL

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


175 kDa

Uniprot ID
[P35573](#)
NCBI
[P35573](#)
Dilution Range

For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.




Western Blot at 1:8000 dilution + human ...

BRAIN
HEART
KIDNEY
LIVER
MUSCLE



Western blot using anti-AGL antibody at ...



Expression of myc-GS causes wild type bu...