


Product Datasheet

AGL Antibody (orb1264728)

Description	AGL Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	This AGL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 357-387 amino acids from the Central region of human AGL.
Target	AGL
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	175 kDa
Uniprot ID	P35573
NCBI	P35573
Dilution Range	For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of hAGL-C371 in K5...



Confocal immunofluorescent analysis of A...