


Product Datasheet

CSRNP1 Antibody (orb1264759)

Description	CSRNP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This CSRNP1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 501-533 amino acids from human CSRNP1.
Target	CSRNP1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	64 kDa
Uniprot ID	Q96S65
NCBI	Q96S65
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western Blot at 1:1000 dilution + human ...