
Product Datasheet

PFN1 Antibody (orb1264780)

Description

PFN1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, WB

Immunogen

This Profilin-1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 108-140 amino acids from the human region of human Profilin-1.

Target

PFN1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25 For WB starting dilution is: 1:2000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


15 kDa

Uniprot ID
[P07737](#)
NCBI
[P07737](#)
Dilution Range


For FACS starting dilution is: 1:25 For WB starting dilution is: 1:2000

Expiration Date


12 months from date of receipt.



Overlay histogram showing HeLa cells sta...



Overlay histogram showing NIH/3T3 cells ...



Western Blot at 1:2000 dilution Lane 1: ...