

---

## Product Datasheet

### ITIH4 Antibody (orb1264805)

**Description**

ITIH4 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

FC, WB

**Immunogen**

This ITIH4 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 784-816 amino acids from human ITIH4.

**Target**

ITIH4

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For FACS starting dilution is: 1:25  
For WB starting dilution is: 1:2000

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


103 kDa

**Uniprot ID**
[Q14624](#)
**NCBI**
[Q14624](#)
**Dilution Range**


For FACS starting dilution is: 1:25  
For WB starting dilution is: 1:2000

**Expiration Date**

12 months from date of receipt.



Overlay histogram showing U-2OS cells st...



Western Blot at 1:2000 dilution + human ...