


Product Datasheet

ACO1 Antibody (orb1264842)

Description	ACO1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This ACO1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 354-388 amino acids from the Central region of human ACO1.
Target	ACO1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	98 kDa
Uniprot ID	P21399
NCBI	P21399
Dilution Range	For WB starting dilution is: 1:2000
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...