

Product Datasheet

UBQLN1 Antibody (orb1264938)

Description

UBQLN1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This Ubiquilin1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 296-326 amino acids from the Central region of human Ubiquilin1.

Target

UBQLN1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:100For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

63 kDa

Uniprot ID

[Q9UMX0](#)

NCBI



[Q9UMX0](#)

Dilution Range


For IF starting dilution is: 1:100For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Immunofluorescent analysis of HeLa cells...



Western blot analysis in SH-SY5Y cell li...



Formalin-fixed and paraffin-embedded hum...