





Product Datasheet

UBQLN2 Antibody (orb1264941)

Description		UBQLN2 Antibody	
Species/Host	Rabbit		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	IHC-P, WB		
Immunogen	This Dsk2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 20-50 amino acids from the N-terminal region of human Dsk2.		
Target	UBQLN2		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		Formalin-fixed and paraffin-embedded hum...
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	66 kDa		Formalin-fixed and paraffin-embedded hum...
Uniprot ID	Q9UHD9		
NCBI	Q9UHD9		
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100		
Expiration Date	12 months from date of receipt.		