


Product Datasheet

SEC11A Antibody (orb1265192)

Description	SEC11A Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This SC11A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 95-123 amino acids from the Central region of human SC11A.	Western Blot at 1:1000 dilution + MCF-7 ...
Target	SEC11A	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Western blot analysis in HepG2 cell line...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	21 kDa	
Uniprot ID	P67812	
NCBI	P67812	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	