

Product Datasheet

PIH1D1 Antibody (orb1265421)

Description	PIH1D1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This PIH1D1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 163-192 amino acids from the Central region of human PIH1D1.	Western Blot at 1:1000 dilution + A375 w...
Target	PIH1D1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	32 kDa	
Uniprot ID	Q9NWS0	
NCBI	Q9NWS0	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	