

Product Datasheet

URI1 Antibody (orb1265432)

Description	URI1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This RMP (URI) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 347-376 amino acids from the C-terminal region of human RMP (URI).
Target	URI1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	O94763
NCBI	O94763
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Antibody is used in Western blot to dete...