


Product Datasheet

CHST2 Antibody (orb1265615)

| | | |
|----------------------------|--|---|
| Description | CHST2 Antibody |  |
| Species/Host | Rabbit | |
| Reactivity | Human, Mouse | |
| Conjugation | Unconjugated | |
| Tested Applications | WB | Western blot analysis in human placenta ... |
| Immunogen | This CHST2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 493-522 amino acids from the C-terminal region of human CHST2. | |
| Target | CHST2 | |
| Preservatives | Supplied in PBS with 0.09% (W/V) sodium azide. | |
| Form/Appearance | Liquid |  |
| Concentration | batch dependent | Western blot analysis in mouse cerebellu... |
| Storage | Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures. | |
| Note | For research use only | |
| Application notes | For WB starting dilution is: 1:1000 | |
| Isotype | Rabbit Ig | |
| Clonality | Polyclonal | |
| MW | 58 kDa | |
| Uniprot ID | Q9Y4C5 | |
| NCBI | Q9Y4C5 | |
| Dilution Range | For WB starting dilution is: 1:1000 | |
| Expiration Date | 12 months from date of receipt. | |