


Product Datasheet

NUBP1 Antibody (orb1265850)

Description	NUBP1 Antibody	 Western blot analysis in 293 cell line l...	
Species/Host	Rabbit		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	WB		
Immunogen	This NUBP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 185-212 amino acids from the Central region of human NUBP1.		
Target	NUBP1		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	35 kDa		
Uniprot ID	P53384		
NCBI	P53384		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		