
Product Datasheet

HOXD8 Antibody (orb1265973)

Description

HOXD8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, WB

Applications
Immunogen

This HOXD8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 243-270 amino acids from the C-terminal region of human HOXD8.

Target

HOXD8

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:10~50
For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


32 kDa

Uniprot ID
[P13378](#)
NCBI
[P13378](#)
Dilution Range


For IF starting dilution is: 1:10~50
For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Fluorescent image of HUVEC cell stained ...



Western blot analysis in T47D cell line ...