


---

## Product Datasheet

### ANXA10 Antibody (orb1266079)

<b>Description</b>	ANXA10 Antibody	
<b>Species/Host</b>	Rabbit	
<b>Reactivity</b>	Human, Mouse	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	WB	
<b>Immunogen</b>	This ANXA10 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 274-300 amino acids from the C-terminal region of human ANXA10.	Western blot analysis in NCI-H292 cell l...
<b>Target</b>	ANXA10	
<b>Preservatives</b>	Supplied in PBS with 0.09% (W/V) sodium azide.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	batch dependent	
<b>Storage</b>	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
<b>Note</b>	For research use only	
<b>Application notes</b>	For WB starting dilution is: 1:1000	
<b>Isotype</b>	Rabbit Ig	
<b>Clonality</b>	Polyclonal	
<b>MW</b>	37 kDa	
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q9UJ72">Q9UJ72</a>	
<b>NCBI</b>	<a href="https://www.ncbi.nlm.nih.gov/nuccore/Q9UJ72">Q9UJ72</a>	
<b>Dilution Range</b>	For WB starting dilution is: 1:1000	
<b>Expiration Date</b>	12 months from date of receipt.	