

Product Datasheet

Caspase 3 (Cleaved-Asp175) antibody (orb126608)

Description Rabbit polyclonal antibody to caspase 3 is also known as apopain whose expression can be detected in...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Internal

Concentration 1mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Clonality Polyclonal


Uniprot ID **P42574**

Dilution Range WB: 1:500-3000, IHC-P: 1:50-100, ELISA: 1:20000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from 2...



Immunohistochemical analysis of formalin...