

Product Datasheet

CHMP5 Antibody (orb1266098)

Description	CHMP5 Antibody	MDA-MB231	
Species/Host	Rabbit		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	WB		
Applications	This CHMP5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 178-204 amino acids from the C-terminal region of human CHMP5.	Western blot analysis in MDA-MB231 cell ...	
Immunogen			
Target			CHMP5
Preservatives			Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	25 kDa		
Uniprot ID	Q9NZZ3		
NCBI	Q9NZZ3		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		