


---

## Product Datasheet

### ZNF416 Antibody (orb1266176)

**Description**


ZNF416 Antibody


**Species/Host**


Rabbit


**Reactivity**


Human


**Conjugation**

Unconjugated


**Tested**

WB


**Applications**

Western blot analysis in 293 cell line l...

**Immunogen**

This ZNF416 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 556-582 amino acids from the C-terminal region of human ZNF416.

**Target**

ZNF416

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

67 kDa

**Uniprot ID**
[Q9BWM5](#)
**NCBI**
[Q9BWM5](#)
**Dilution Range**

For WB starting dilution is: 1:1000

**Expiration Date**

12 months from date of receipt.