

Product Datasheet

PVRL4 Antibody (orb1266536)

Description	PVRL4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This PVRL4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 239-266 amino acids from the Central region of human PVRL4.
Target	PVRL4
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q96NY8
NCBI	Q96NY8
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis in A2058 cell line...