


Product Datasheet

ARHGEF18 Antibody (orb1266543)

Description	ARHGEF18 Antibody	 Western blot analysis in 293 cell line l...
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This ARHGEF18 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 91-119 amino acids from the N-terminal region of human ARHGEF18.	
Target	ARHGEF18	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	131 kDa	
Uniprot ID	Q6ZSZ5	
NCBI	Q6ZSZ5	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	